

ICICLE

Thermal Electric Ligature Director

**MELTS
THERMAL
WIRE!**



Prevents Bracket Pop-Offs!

**Nothing You'll Try Will
Work As Well As The Icicle...
*Guaranteed!***

Made By:

ALASKAN
ORTHODONTICS®

- **ALWAYS READY** - One flick of the switch and in 90 seconds it's ready.
- **SAVES PREPARATION TIME** - No filling, spraying or freezing.
- **SAVES CHAIR TIME** - Less time ligating, minimizes bracket pop-off.
- **SAVES MONEY** - Use larger size wires earlier.



***Maintains a Constant
Cold Temperature
of 36°***

Within seconds of applying the tip to the archwire, the wire becomes soft and limp, enabling the director slot-shaped tip to easily lift the wire into the bracket slot.

EASY STERILIZATION

1. Remove the tip using an instrument. replace with a clean one, dry heat or autoclave.
2. Remove the used cover, dispose of and replace with a clean one.

Portable unit shown also available built into most Ross Equipment.

The Icicle measures 6"x 6-1/4"x 2-1/2 and requires a 30-psi air supply and a 110-volt outlet.

ROSS

ORTHODONTIC

INDEX