

# Flexion Stool

*Ergonomically  
Designed*

*Posturpedic  
Seat*

*Double  
Contoured*



The Flexion Stool incorporates a molded plywood seat pan and backrest with compound contours that help keep the stool with you as you move – *unlike round or flat stools you have to fight to stay on!*

**ROSS**  
ORTHODONTIC

MAXUS CHAIR | OPTIMA CHAIR | OPTIMA CHAIR WITH POWER BASE

FLEXION ASSISTANT'S STOOL | SIMPLICITY STOOL | INDEX

# FLEXION STOOL

## A Good Stool Will Take Unnecessary Forces Off Your Body

*A good stool unloads the spine and keeps your body in a "neutral" position – your ear, shoulder and hip are in vertical alignment when your feet are on the floor.*



**Resigning yourself to settle for less is like settling for a shoe that is a half-size too small; you don't always feel it until the end of the day.**

### ***Double Contoured Seat***

The steep waterfall of the seat promotes good leg circulation. The seat pan is fixed at a -2° angle for comfort when the operator is seated for a long period of time.

### ***Backrest Adjustment***

With a twist of a knob, you can adjust the backrest height. The backrest horizontal location is spring-loaded and can be adjusted by a pull of a lever and locked in position for good lumbar support.

### ***Automatic Height Adjustment***

The pneumatic height adjustment, activated by a soft pull lever, adjusts from 17" to 22".

### ***Stability and Strength***

The 5-legged, one-piece cast aluminum base is reinforced with cast-in welts for extra stability and easy movement.

### ***Thick Seat and Back Padding***

The seat and back padding are made from 2" of foam, 1" of firm foam and 1" of medium foam, for assured comfort when seated for long periods of time.

### ***Available in over 150 Standard Colors***

Boltaflex Color Guard and Grand Sierra with Prefix®

### ***Optional: Ultra Leather***

**ROSS**

ORTHODONTIC

1-800-247-4109 • 1-972-775-8757 • FAX 1-972-775-8758 • E-Mail rossortho@att.net

PO Box 882, Midlothian, TX 76065-0882

**Office Design and Task Orientated Equipment**