

Protegé

This unit can be located on the right side of the chair as shown, or a mirror version is available for the left side.



Unit: 23" W X 27" D X 31" H
Drawers: 8 1/2" W X 12" D X 2 3/4" H
Unit With Ross Chair 48"

**Chairside Storage • Large Work Surface
Trans-Thorax or Side Delivery • Radius Corners
Aseptic Design**

ROSS
ORTHODONTIC

ROSS MODEL A | SATURN II MOBILE | MINI PRO MOBILE | INDEX

PROTEGE

Instrumentation

High volume evacuation assembly, saliva ejector assembly, solids collector, three-way quick-change syringe assembly, two hand piece control, one high speed (wet), one slow speed (dry), disc type foot control, air filter/regulator, master on-off switch, and bottled water system.

Chairside work area

The unit top offers over 470 square inches of chairside work area.

Chairside storage

Durable drawers constructed of Baltic Birch, UV sealed and installed using 100 lbs. full extension ball-bearing glides for smooth operation and full access. Over 2,140 cubic inches of chairside storage space.

Durability

The unit is constructed of 3/4" cabinet-grade plywood for durability and retention of integrity for the life of a practice. The steel sub-base ensures a solid mounting to the floor for long-term stability.

Options

- Clean water system
- Computer monitor shelf
- Pull out work tablet
 - Chart holder



O R T H O D O N T I C

1-800-247-4109 • 1-972-775-8757 • FAX 1-972-775-8758 • E-Mail rossortho@att.net

PO Box 882, Midlothian, TX 76065-0882

Office Design and Task Orientated Equipment